


Some thoughts on the benefits of including Occlusograms in the Cephalometric analysis of Growth and Treatment.

The original decision to include Occlusograms in the Tiops program was inspired by the ideas of Professor Björk in connection with 'Cephalometric Analysis of Growth and Treatment', and are described in: <http://tiops.com/downloads/Articles/BjorkCephalometricGrowthAnalysis.pdf> from where this illustration originates:



What can we learn by including the occlusogram?


1. One of the most important benefits of including occlusograms of the study casts is that we now can study the occlusal changes in all three dimensions. Information about the sagittal (horizontal) and vertical changes of the dentition, we can get from the lateral headfilm, but the transverse changes only from the occlusogram.
2. By including the occlusograms in the cephalometric analysis, we can get a more precise idea of the transverse changes of both the posterior and anterior teeth. Expansion of the posterior teeth and midline changes of the anterior teeth can now also be more precisely determined.



About midlines


To what degree the clinician chooses to incorporate correction of the midline in the treatment plan is a personal choice. However, it is important to analyze the midline situation in order to be able to make this decision. An important part of the analysis is to distinguish between the Facial and Dental midlines, as seen in the Frontal Aspect and the Dental Arch Midline or Symmetry Axis.

The patient's Facial Midline, as seen in the Frontal Aspect, is best diagnosed using a Frontal Photo or a Frontal Cephalometric Headfilm. It is important that the face is correctly aligned without any rotation.




About Arch Analysis Midline Corrections


In relation to performing an Arch Analysis it is usually a good idea to check the **Arch Midline or Symmetry Axis** placement before making any corrections to the molar position (as any "wrong" molar position may be due to an incorrectly located midline). Do not spend too much energy on precisely defining the midline at the beginning of the registration procedure as it is easier to correct the midline after all teeth have been digitized. To make a necessary correction of the midline placement, use the function **View/Show Symmetry Guides**, and open the Points list under Upper/Lower Arch and scroll down to the landmark ref1




Before correction



During correction



After correction





Example showing the incisal midline not coinciding with the symmetry axis

About correcting the molar position on the cephalometric headfilm using the dental arch tracings

If the initially digitized molar position on the Lateral Analysis does not match the position on the occlusogram, it can be corrected using the following procedure.

When selecting the function **View/Show X-Ray Teeth on Arches**, guide Lines showing the mesial aspects of the Upper and Lower **Arch Analysis Molars** will now automatically appear in the **Lateral Analysis** displayed in the Main Window. This is in addition to the shaded Lateral Analysis Molars, that appear in the Jaw View Window. This function works in combination with the function **View/Show Right/Left Side** that make separate corrections of either side possible. This function works when the molar on only one side has been registered. However, it also works when both the Right and Left Molars are included in the Lateral Analysis.



Before correction:
The Cephalometric Upper Right Molar is placed too far distally in relation to the Right Molar at the Arch Analysis

In order to make the necessary correction open the **Points** list under Xray and scroll down to the landmark **msm**.

Stage1	
[- Xray]	
Filename	C:\TIOPS\Analy4\Xray\
Mirror	False
ResolutionX	203
ResolutionY	203
Enlargement	5.6
[- Filters]	
Invert	False
Brightness	1
Contrast	-53
SharpenPercent	0
[- Photo]	
[- Points]	
re1	[-9489,15049]
re2	[-9114,-13385]

isi	[-4631,-7626]
asi	[-4617,-5198]
msm	[-2657,-6701]
rsm	[-2560,-5274]

To maneuver the tooth:
 A. MouseMove to rotate
 B. Shift+MouseMove to translate


Procedure:

1. Place the tooth correctly using A and B
2. LeftClick to lock in position


The tooth will always be registered in the equal size irrespective of the cursor position along the tooth axis

msm
rsm

Tooth Positioning Procedure
 Upper Molar



When the molar is corrected in the **Main Window**, it is automatically updated in the **Jaw View** window.



About Printouts



The Tiops Program offers four different printout options of the Arch Analysis.

1. **Upper or Lower Arch** with one or more Stages.
2. **Upper/Lower Arch – OneStage** sheet.
3. Occlusograms superimposed on different reference landmarks called the **Upper/Lower Arch – MoreStages – 2 Superimpositions** sheet.
4. Printout of the Arch Variables

To print the appropriate sheet(s) just select the desired stages and open the actual analysis in the Main window.

Note: Any Enlargement factor can be selected.

1. Printout of the **Upper or Lower Arch** including one or more Stages.




Enlargement: 0%

Stage 1 (MA)



- 07-12-2007 SkelAge 14:3




Stage 2 (TA)

- 09-04-2010 SkelAge 16:8



2. Printout of **Upper/Lower Arch – OneStage**





Note: Upper and Lower Arch superimposition using the cephalometric relations

r6 MesDis	mm	10.0	10.3	-0.4
r5 MesDis	mm	7.1	7.1	-0.1
r4 MesDis	mm	6.9	6.8	0.1
r3 MesDis	mm	6.2	6.6	-1.5
r2 MesDis	mm	5.6	5.9	-0.9
r1 MesDis	mm	5.4	5.5	-0.3
l1 MesDis	mm	5.4	5.5	-0.4
l2 MesDis	mm	5.5	5.9	-1.2
l3 MesDis	mm	6.0	6.6	-2.0
l4 MesDis	mm	6.5	6.6	-0.5
l5 MesDis	mm	6.5	7.1	-1.3
l6 MesDis	mm	10.3	10.3	0.1


r6 MesDis	mm	9.2	10.2	-1.7
r5 MesDis	mm	6.2	6.5	-0.6
r4 MesDis	mm	6.8	6.6	0.4
r3 MesDis	mm	7.5	7.5	-0.1
r2 MesDis	mm	7.0	6.8	0.4
r1 MesDis	mm	8.3	8.7	-0.6
l1 MesDis	mm	8.3	8.7	-0.6
l2 MesDis	mm	7.0	6.8	0.3
l4 MesDis	mm	6.6	6.6	0.0
l5 MesDis	mm	6.4	6.5	-0.1
l6 MesDis	mm	9.7	10.2	-0.8

LoAntDepth	mm	9.0	13.0	-3.1
LoAntWidth	mm	32.5	31.5	0.6
LoPostWidth	mm	50.0	44.8	2.1
Space1To5R	mm	-1.2	0.0	-0.6
Space1To5L	mm	0.1	0.0	0.1
Pred Space1To5R	mm	-1.6	0.0	-0.8
Pred Space1To5L	mm	-1.4	0.0	-0.7

LoAntDepthDiscr	mm	2.2	0.0	1.2
LoAntWidthDiscr	mm	2.7	0.0	0.9
LoPostWidthDiscr	mm	6.1	0.0	2.0

UpAntDepth	mm	9.8	16.0	-4.5
UpAntWidth	mm	29.8	31.5	-1.2
UpPostWidth	mm	44.0	44.8	-0.3
Space1To5R	mm	-3.9	0.0	-1.9
PredSpace1To5R	mm	-4.4	0.0	-2.2
PredSpace1To5L	mm	-4.5	0.0	-2.2

3. Printout of occlusograms superimposed on different reference landmarks called the **Upper/Lower Arch – MoreStages – 2 Superimpositions** sheet.



Print - LowerArch [X]

General

Select Printer
Printer: Canon TS9100 series

Sheets

Tracing Structured Variable
 Variables Cephalometric Coordinates
 Explanation Full Head Tracing

Enlargement: 0

Lateral - MoreStages - NoArches - Bold
 Upper/Lower Arch - MoreStages - 2 Superimpositions

PrinterSetup Print Cancel Help

Enlargement: 5.6%



Stage 1 (MA)

- 07-12-2007 SkelAge 14:3

Stage 2 (TA)

- 09-04-2010 SkelAge 16:8

Overlay at Model Incisor



Overlay at Cephalometric Incisor

