

Stage types

Type	Description	Automatically generated Landmarks	
		Landmarks	Reference Points
MA	Morphological Analysis		ma1 - ma2 - mx1 - mx2
GP	Growth Progress	spg	ma1 - ma2 - mx1 - mx2
TP	Treatment Progress	spg	ma1 - ma2 - mx1 - mx2
TA	Treatment Analysis	spg	ma1 - ma2 - mx1 - mx2
RP	Retention Progress	spg	ma1 - ma2 - mx1 - mx2
SA	Stability Analysis	spg	ma1 - ma2 - mx1 - mx2
MS	Morphological Surgical Analysis		ma1 - ma2 - mx1 - mx2 - nl - nl1 - nl2
BS	Before Surgery	spg	ma1 - ma2 - mx1 - mx2 - nl - nl1 - nl2
PS	Progress Surgery	spg - pm - pal - ss - sp - teeth	ma1 - ma2 - mx1 - mx2 - nl - nl1 - nl2
SA	After Surgery Stage	spg - pm - pal - ss - sp - teeth	ma1 - ma2 - mx1 - mx2 - nl - nl1 - nl2
SS	Surgical Stability Analysis	spg	ma1 - ma2 - mx1 - mx2 - nl - nl1 - nl2